PolyMem® WIC®-sårfyldstoff

BESKRIVELSE:

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PolyMem® WIC® Wound Filler

DESCRIPTION:

PolyMem WIC wound filler products are primary dressings which may be placed into open wounds to eliminate dead space, absorb exudate, and help establish and maintain the appropriate moisture balance of exposed tissues. They provide fast wicking, gently expand (about 1/3) to fill space and help create an optimal healing environment. The filler minimizes the need to disturb the wound bed. It allows extended time between dressing changes and provides easier changes. Fewer and faster changes save time and money. PolyMem WIC products are made of moisture loving polyurethane matrix that contains (1) safe wound cleanser, (2) glycerin moisturizer, and (3) superabsorbents.

PolyMem WIC Silver products incorporate very fine silver particles, and will not stain skin. Silver acts inside the dressing to protect from microbial contamination. Approximate silver content by area 186 micrograms/cm² in the formulation.

Drug-free PolyMem WIC helps to:

- effectively manage and heal wounds
- absorb fluid and provide a moist healing environment
- minimize and relieve persistent and procedural wound pain
- reduce edema, bruising, and the spread of inflammation into surrounding undamaged tissues

PolyMem WIC Cavity Fillers are perforated in strips to facilitate ease of use (for separation or folding). PolyMem WIC Silver Rope may be used in tunneling wounds. The top and bottom wide surfaces of the dressing are reinforced so that the entire dressing may be easily removed when placed in tunnels. All PolyMem dressings help relieve wound pain.

INDICATIONS FOR ALL HUMAN BEINGS:

Under the direction of a healthcare professional, PolyMem wound fillers may be used for the management of:

- Pressure ulcers (stages III-IV)
- Diabetic Ulcers
- Vascular Ulcers
- Acute Wounds

Any serious incident that has occurred in relation to this device should be reported, by the user and/or patient, to Feris Mfg. Corp. and the competent authority of the Member State in which the user and/or patient is established.

PRECAUTIONS:

- Use only in wounds where there is sufficient room to contain the wound filler and to facilitate removal later. Trim so the filler is about 1/3 smaller than the open cavity or tunnel width. Leave rope longer than the tunnel to make it easier to remove.
- Do not pack tightly
- For external use only
- PolyMem dressings are not compatible with oxidizing agents (hydrogen peroxide or hypochlorite solutions)
- Prior to performing imaging procedures, determine if the dressing should be removed.
- Avoid contact with electrodes or conductive gels during electronic measurements (e.g. EEG or ECG).
- Topical treatments are not recommended in conjunction with PolyMem dressings. Use of PolyMem dressings with lotions, creams, or ointments may diminish the benefits of the PolyMem dressing or change the effect of the other products.
- Do not use and discontinue use on people who show signs of sensitivity, irritation, or allergy from the dressing or its materials
- Intended for one use on one person. Do not reuse. Reuse may degrade properties, spread contamination, and increase risk of infection.

PolyMem dressings may be used when signs of infection are present if proper medical treatment that addresses the cause of the infection has been implemented.

PolyMem WIC Silver dressings are recommended when a wound displays delayed healing, signs of infection, or increased bioburden; or the patient may display infection risk factors (such as contaminated wound, poor nutrition, diabetes, history of infection, etc.)

Store at room temperature ≤30°C (86°F).

Excursions permitted <55°C (131°F) but should be minimized. Dispose as typical for biohazardous waste.

Dimensions are approximate.

DIRECTIONS:

Eliminate cavities by gently filling, AVOID OVERFILLING THE WOUND. Over-packing may increase pressure on the wound, potentially causing damage. Use sterile or clean technique as appropriate for the situation.

APPLICATION:

1. Initial dressing application: Prepare the wound according to protocol or as directed by a physician or ordering clinician.
2. Size - Cavity filler and rope will expand about 1/3 when absorbing fluid. Therefore, they should be smaller and thinner than the wound by 1/3. For open wounds, Cavity Filler is perforated in 1" wide strips for easy folding or detachment, or may be cut to size. For tunneling wounds, Rope should be narrower than the tunnel or wound.


PolyMem WIC can be divided along perforations

Perforations can also aid in folding WIC

If narrower Rope is desired, cut in half using the slits at each end as a guide.

4. Cover - PolyMem WIC products are primary dressings. Cover with suitable PolyMem dressing as a secondary dressing. Helpful hint - Draw the wound size on the outside of the secondary dressing to help determine when a change is needed.

Dressing Change:

1. For an exuding wound, fluid will become visible through the top of the PolyMem secondary dressing. Change before fluid reaches the wound margin, when good practice dictates, or after no more than 7 days. Change immediately if fluid reaches the edge of the covering dressing.

2. As with other dressings, more frequent changes may be indicated if the patient has a compromised immune system, diabetes, infection or when dressing to speed up the removal of non-viable tissue in the wound.

3. To remove, lift and pull slowly and gently by hand or with wide tip forceps. Small, sharp forceps may cut or shred saturated filler. PolyMem WIC fillers will not adhere to the wound bed, usually assuring pain-free changes. However, inspect the dressing and wound to assure all material has been removed.

4. If coagulated blood or dried exudate has caused the dressing to adhere to skin or hair, use a small amount of saline or water to slowly gently soften and loosen. Clean the intact skin around the wound.

5. Manage dressing changes according to protocol or as directed by a physician or other ordering clinician. In most cases, when using PolyMem, there is no need to disturb or cleanse a wound during changes unless the wound is infected or contaminated. PolyMem is designed to continuously clean the wound and does not leave residue that needs to be removed.

Excessive cleaning may injure regenerating tissue and delay wound healing.

Observations during use or change:

- A dramatic increase in wound fluid may be observed during the first few days due to the fluid attraction. This is not uncommon and indicates that the dressing is working.
- During normal healing, the wound may appear to become larger during the first few dressing changes because the body is naturally removing non-viable tissue. If the wound continues to become larger, consult an expert.
- PolyMem WIC dressings help support new blood vessel formation so it is not uncommon to see blood-saturated wound fluid and dressings during changes.
- Be alert for signs and symptoms, such as, increased redness, inflammation, pain, warmth, odor, whitening, fragility, excessive